

Manufacturer	Implant	Torque (Ncm)	Torque Settings	Driver & Head
BioHorizons	Internal 3.5, 4.5, 5.7	30 Ncm	0.050" (1.25 mm) Hex	⬢
	Tapered Internal 3.5, 4.5, 5.7	30 Ncm	0.050" (1.25 mm) Hex	⬢
	External Implant 3.5	30 Ncm	0.050" (1.25 mm) Hex	⬢
	External Implant 4.0, 5.0	30 Ncm	0.050" (1.25 mm) Hex	⬢
Dentsply	ANKYLOS® C/X – C/ 3.5, 4.5, 5.5, 7.0 (non-indexed)	15 Ncm	1.0 mm Hex	⬢
	ANKYLOS® C/X – /X 3.5, 4.5, 5.5, 7.0 (indexed)	15 Ncm	1.0 mm Hex	⬢
	OsseoSpeed™ EV 3.0 – Green	25 Ncm	0.050" Hex	⬢
	OsseoSpeed™ EV 3.6 – Purple	25 Ncm	0.050" Hex	⬢
	OsseoSpeed™ EV 4.2 – Yellow	25 Ncm	0.050" Hex	⬢
	OsseoSpeed™ EV 4.8 – Blue	25 Ncm	0.050" Hex	⬢
	OsseoSpeed™ EV 5.4 – Brown	25 Ncm	0.050" Hex	⬢
	OsseoSpeed™ Profile EV 4.2 Yellow	25 Ncm	0.050" Hex	⬢
	OsseoSpeed™ Profile EV 4.8 Blue	25 Ncm	0.050" Hex	⬢
	OsseoSpeed™ / OsseoSpeed™ TX 3.0 – Yellow	15 Ncm	0.050" Hex	⬢
	OsseoSpeed™ / OsseoSpeed™ TX 3.5, 4.0 – Aqua	20 Ncm	0.050" Hex	⬢
	OsseoSpeed™ / OsseoSpeed™ TX 4.5, 5.0 – Lilac	25 Ncm	0.050" Hex	⬢
	OsseoSpeed™ 4.0 - Blue (Previous)	25 Ncm	0.050" Hex	⬢
	OsseoSpeed™ TX Profile 4.5, 5.0	25 Ncm	0.050" Hex	⬢
	XIVE® S 3.0	24 Ncm	1.22 mm Hex	⬢
	XIVE® S 3.4, 3.8, 4.5, 5.5	24 Ncm	1.22 mm Hex	⬢
FRIALIT™ 3.4, 3.8, 4.5, 5.5	24 Ncm	1.22 mm Hex	⬢	
Nobel Biocare®	Angulated Screw Channel	35 Ncm	Nobel Biocare® Omnigrip	⬢
	NobelActive™ 3.0	15 Ncm	Nobel Biocare® Unigrip	★
	NobelActive™/NobelReplace™ Conical Connection NP 3.5	35 Ncm	Nobel Biocare® Unigrip	★
	NobelActive™/NobelReplace™ Conical Connection RP 4.3, 5.0	35 Ncm	Nobel Biocare® Unigrip	★
	NobelReplace™ NP 3.5	35 Ncm	Nobel Biocare® Unigrip	★
	NobelReplace™ RP 4.3	35 Ncm	Nobel Biocare® Unigrip	★




Implant Torque Specification Guide

Ensure you are applying the right amount of torque for your patient's implant restoration courtesy of Burdette Dental Lab.


Questions? Call Burdette at 1-800-624-5301

Impression Techniques - Closed Tray Instructions


Step 1
Remove healing abutment




Step 2
Insert closed-tray coping—OR—
Use plastic impression tray



Step 3
Confirm proper seating




Step 4
Block out driver indentation




Step 5
Inject impression material



Step 6
Replace healing abutment



Step 7
Place replica in impression




Manufacturer	Implant	Torque (Ncm)	Torque Settings	Driver & Head
Nobel Biocare®	NobelReplace™ WP 5.0	35 Ncm	Nobel Biocare® Unigrip	★
	NobelReplace™ 6.0	35 Ncm	Nobel Biocare® Unigrip	★
	Brånemark System™ NP 3.3	35 Ncm	Nobel Biocare® Unigrip	★
	Brånemark System™ RP 3.75, 4.0	35 Ncm	Nobel Biocare® Unigrip	★
	Brånemark System™ WP 5.0, 5.5	35 Ncm	Nobel Biocare® Unigrip	★
Straumann®	Bone Level 3.3 NC (Narrow CrossFit™)	35 Ncm	SCS Screwdriver	★
	Bone Level 4.1, 4.8 RC (Regular CrossFit™)	35 Ncm	SCS Screwdriver	★
	Standard/Standard Plus 4.8 RN (Regular Neck)	35 Ncm	SCS Screwdriver	★
	Standard/Standard Plus 4.8 WN (Wide Neck)	35 Ncm	SCS Screwdriver	★
	Standard Plus 3.3 NN (Narrow Neck)	35 Ncm	SCS Screwdriver	★
ZimVie Dental	Tapered Screw-Vent 3.5, 4.5, 5.7/ Tapered Screw-Vent 3.5, 4.5, 5.7 friction fit	30 Ncm	1.25 mm Hex	⬢
	Screw-Vent 3.3, 3.7, 4.5/ Screw-Vent 3.3, 3.7, 4.5 friction fit	35 Ncm	1.25 mm Hex	⬢
	SwissPlus™ 4.8	35 Ncm	1.25 mm Hex	⬢
	Certain™ 3.25, 4.0, 4/3, 5.0, 5/4, 6.0	20 Ncm	0.048" Hex	⬢
	Certain™ Prevail™ 3/4/3 (3.4), 4/5/4, 5/6/5, 4/3, 5/4, 6/5	20 Ncm	0.048" Hex	⬢
	Certain™ XP 4/5, 5/6	20 Ncm	0.048" Hex	⬢
	Miniplant™ 3.25	20 Ncm	0.048" Hex	⬢
	Osseotite™ 3.75, 4.0, 5.0, 6.0	20 Ncm	0.048" Hex	⬢
	Osseotite™ XP 4/5, 5/6	20 Ncm	0.048" Hex	⬢
Megagen	MegaGen	35 Ncm	0.048" (1.2 mm) Hex	⬢
	MegaGen Anyridge	35 Ncm	0.048" (1.2 mm) Hex	⬢
	Mini	15 Ncm	0.048" (1.2 mm) Hex	⬢

Please note: When using multi-unit abutments, it is recommended they be torqued to no more than 15 Ncm.

Implant Tool Loaner Program

It's difficult for any practice to maintain all the tools necessary to work on every platform. With our Implant Tool Loaner Program, now doctors are no longer limited in what they can prescribe and work on. Burdette can provide the exact equipment you need to seat our abutments on any implant—at no additional cost!



Questions? Call Burdette at 1-800-624-5301

Impression Techniques - Open Tray Instructions

