



131 Lyon Lane • Birmingham, AL 35211 • 800-624-5301 • 205-916-0887
info@burdettedental.com • www.burdettedental.com

Doctor Name: _____ Phone Number: _____

Email: _____

Patient Name: _____ Patient Age: _____ Gender: M F

Deliver case before appointment.

Deliver case based on lab capacity.

Date: _____ by 5 P.M.

(Patient will be appointed after receiving case)

Instructions

\$79 Advanced Z16™ and \$89 Advanced Z16™ ML

*Price per unit up to three units.

\$79 Per Unit
on Advanced Z16™

\$89 Per Unit
on Advanced Z16™ ML

Need Supplies?

- Boxes
- Bio Bags
- FedEx Shipping Labels
- RX Forms
- Other: _____

Enclosed with Case

- Impression Triple Tray
- Study Model Opposing
- Bite
- Photos
- Other: _____

Shade Instructions

Shade: _____



Signature: _____

License #: _____ Date: _____ Promo Code: DPRADVQ121

Payment is due upon receipt of statement. Total statement amount due by the end of the month. All past due invoices will be subject to a finance charge and/or collection fees. The undersigned is responsible both corporately and personally. Your signature is acceptance of these terms.

For Lab Use Only

Contact/LMOM: _____ Date: _____


RE: _____

Result: _____

Initial: _____

Preparation Guidelines

Ideally, the preparation includes a circumferential shoulder or chamfer with a horizontal angle of at least 5° . The vertical preparation angle should be at least 4° . The inside angle of the shoulder preparation must be given a rounded contour. All occlusal and incisal edges should also be rounded. The marginal edge of the preparation needs to be continuous and clearly visible. A bevel should be avoided. For posterior and anterior teeth, a supragingival margin poses no problems. Due to the tooth-colored framework, very esthetic results can be achieved.



Do:

- Change nothing from your usual crown and bridge/PFM technique. Prep, imp, bite, cement.
- As with any ceramic restoration, proper reduction of tooth/prep helps us build a restoration of unequal esthetics.

Don't:

- Leave sharp tips on corners of the prep.
- Do a shoulder prep. This is not necessary for Crystal or Diamond crowns. If you do prep a shoulder, please do not prep a trough on the shoulder (the scanner does not read these well).
- Have patient force crown to place with cotton roll or anything else.
- Internally relieve all ceramic crown to seat. This will cause micro-fractures and in return will cause breakage in the future. If you feel relief needs to be done for additional seating, relief must be done to prep.

Terms & Policies

Payment is due upon receipt of statement. Total statement amount due by the end of the month. All past due invoices will be subject to a finance charge and/or collection fees. The undersigned is responsible both corporately and personally. Your signature is acceptance of these terms.