



iTero® Element

Connect Your iTero® to Burdette's Workflow

Save time and money when you send scans directly from your iTero® digital impression scanner. **In addition to all-ceramic restorations, Burdette accepts digital impressions for PFM and full-cast restorations.**

- Call iTero® tech support at 800-577-8767, select option 2 and ask them to add Burdette as a lab partner.
- Select Burdette as your lab partner and submit your digital file.

If you have any questions, contact Burdette at 800-624-5301 and ask for digital impression technical support.

Looking for More Information on Our Products?

For material descriptions, preparation and cementation guides and billing codes, visit burdettedental.com/products.

Save money on every digital submission by choosing model free when possible.

SAVE \$20 ON YOUR FIRST DIGITAL CASE

Digital submissions are simple when you choose Burdette. Reference **promo code: GO_DIGITAL** on your Rx when you upload your first case and save \$20 on your first digital submission.

Promo Code: GO_DIGITAL