



## Connect Your iTero® to Burdette's Workflow

Save time and money when you send scans directly from your iTero® digital impression scanner. **In addition to all-ceramic restorations, Burdette accepts digital impressions for PFM and full-cast restorations.**

- Call iTero tech support at 800-577-8767, select option 2 and ask them to add Burdette as a lab partner.
- Select Burdette as your lab partner and submit your digital file.

If you have any questions, contact Burdette at 800-624-5301 and ask for digital impression technical support.

### Looking for more information on our products?

For material descriptions, preparation and cementation guides and billing codes, visit [burdettedental.com/products](http://burdettedental.com/products)

Save money on every digital submission by choosing model free when possible.

### SAVE \$20 On Your First Digital Case

Digital submissions are simple when you choose Burdette. Reference coupon code **GO\_DIGITAL** on your Rx when you upload your first case and save \$20 on your first digital submission.

**Coupon Code: GO\_DIGITAL**